English

Deutsch

Français

Italiano

Español


Operating Instructions **FiveEasyPlus™** pH meter FEP20

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METTLER **TOLEDO**

Operating Instructions
FiveEasyPlus™

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Español

Table of contents

1	Introduction	3
2	Safety measures	4
3	Installation	6
3.1	Unpacking	6
3.2	Installing the Electrode Arm (optional)	6
4	Display and Key controls	7
5	Operating the Instrument	9
5.1	Calibration	9
5.1.1	Buffer Groups	9
5.1.2	Selecting a predefined buffer group	9
5.1.3	Performing a 1 or 2-Point Calibration	9
5.1.4	Performing a 3-Point Calibration	9
5.2	Sample Measurement	10
5.2.1	Performing a pH Measurement	10
5.2.2	Performing a mV / Redox Measurement	10
5.3	Temperature Measurement	10
5.3.1	Manual Temperature Compensation	11
5.4	Using the Memory	11
5.4.1	Storing a Reading	11
5.4.2	Recalling from Memory	11
5.4.3	Clearing the Memory	11
5.5	Printing	11
5.5.1	Connection and Configuration	11
5.5.2	Printing a Measurement Results	12
5.5.3	Printing Calibration Results	12
5.5.4	Printing from Memory	13
5.6	Self-Diagnosis	13
5.7	Error Messages	13
6	Maintenance	15
6.1	Meter Maintenance	15
6.2	Electrode Maintenance	15
6.3	Disposal	15
7	Accessories	16
8	Specifications	17
9	Appendix	18
9.1	Buffer Tables	18
9.2	Error Limits	19

1 Introduction

Thank you for purchasing this high quality METTLER TOLEDO meter. Ease of use, reliable results and ergonomic design – these are our goals for products.

The FiveEasy meters have an excellent price/performance ratio and many useful accessories are already included in the standard delivery of the instruments. Other accessories such as electrodes, solutions or electrode arms are of course also available from METTLER TOLEDO.

You will find much helpful functionality in these meters. A few of many exiting features are:

- Automatic endpoint format which freezes the measurement automatically when the value is stable.
- Four predefined buffer groups to select from.
- Automatic buffer recognition which allows you to calibrate with your buffers in any order.

2 Safety measures

Measures for your protection

Risk of explosion

- Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight (explosion hazard due to spark formation, corrosion caused by the ingress of gases).

Risk of corrosion

- When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules!

Caution

- Never unscrew the two halves of the housing!
- Have the meter serviced only by METTLER TOLEDO Service!
- Any spillage should be wiped off immediately! Some solvents might cause corrosion of the housing.
- Avoid the following environmental influences:
 - Powerful vibrations
 - Direct sunlight
 - Atmospheric humidity greater than 80%
 - Corrosive gas atmosphere
 - Temperatures below 5 °C and above 40 °C
 - Powerful electric or magnetic fields

FCC Rules

This device complies with Part 15 of the FCC Rules and Radio Interference Requirements of the Canadian Department of Communications. Operation is subject to the following conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.


This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

3 Installation

3.1 Unpacking


Carefully unpack the meter. Keep the calibration certificate in a safe place.

Insert the right adapter clip into the power adapter slot:


3.2 Installing the Electrode Arm (optional)

The electrode arm can be fastened to FiveEasyPlus™ on the left or the right side.

- 1 Remove the electrode stand rubber cover on that side of the bench meter where you prefer the stand. Keep the cover in a safe place.
- 2 Insert the electrode stand firmly into the opening as shown.
- 3 Fix the opening of the arm upon the stand from above, adjust it to a moderate height.

4 Display and Key controls

- 1** Electrode condition (For electrode maintenance see "Electrode Maintenance (page 15)")

Slope: 95-105%
and offset: \pm (0-15) mV
Electrode is in good condition

Slope: 90-94%
or offset: \pm (15-35) mV
Electrode needs cleaning

Slope: 85-89%
or offset: \pm (>35) mV
Electrode is defective

- 2** Calibration icon - calibration in progress

- 3** Measurement icon - measurement or calibration is running

- 4** Setup icon - meter is in setup mode

- 5** pH/mV reading or slope in calibration process

- 6** Auto / Manual temperature compensation (see "Temperature Measurement (page 10)")

- 7** Endpoint stability \checkmark , Auto endpoint $/\bar{A}$

- 8** Temperature during measurement or offset in calibration process

- 9** Calibration point **Cal** / Buffer group number **MR** / Memory number **MR** / Error index **Err**

	Press & release	Press & hold for 3 seconds
	<ul style="list-style-type: none"> Start or endpoint measurement Confirm setting, store entered value 	<ul style="list-style-type: none"> Turn auto endpoint on /A / off /
	<ul style="list-style-type: none"> Start calibration 	<ul style="list-style-type: none"> Review the latest calibration data
	<ul style="list-style-type: none"> Meter on Back to measurement screen 	<ul style="list-style-type: none"> Meter off
	<ul style="list-style-type: none"> Store current reading to memory Increase value during setting Scroll up through the memory 	<ul style="list-style-type: none"> Recall stored data Print current data set of stored data
	<ul style="list-style-type: none"> Switch between pH and mV measuring modes Decrease value during setting Scroll down through the memory 	<ul style="list-style-type: none"> Enter Setup mode
	<ul style="list-style-type: none"> Start self-diagnosis 	

5 Operating the Instrument

5.1 Calibration

5.1.1 Buffer Groups

The FiveEasyPlus™ pH meter allow you to perform 1, 2, or 3 point calibrations. After selecting your calibration buffer group from the predefined buffer groups in the meter, the buffers are automatically recognized and displayed during calibration (auto buffer recognition).

The four predefined groups are:

B1	1.68	4.01	7.00	10.01		(at 25°C)
B2	2.00	4.01	7.00	9.21	11.00	(at 25°C)
B3	1.68	4.00	6.86	9.18	12.46	(at 25°C)
B4	1.68	4.01	6.86	9.18		(at 25°C)

Tables (B1...B4) for automatic temperature compensation are programmed in the meter for each buffer group (see also "Appendix (page 18)").

5.1.2 Selecting a predefined buffer group

- 1 Press and hold the **Mode/Setup** key until the setup icon appears on the display and the **MTC** temperature blinks.
- 2 Press **Read** to ignore – the current temperature unit blinks – press **Read** to ignore again.
⇒ When the current buffer group blinks, the buffers belonging to this buffer group appear alternating on the display.
- 3 Use **▲** or **▼** to select another buffer group and press **Read** to confirm your selection.
⇒ The meter will exit automatically to the measurement screen.

5.1.3 Performing a 1 or 2-Point Calibration

- Place the electrode in a calibration buffer and press **Cal**.
⇒ The calibration icon and the measurement icon appear on the display.
The meter endpoints according to the preselected endpoint mode after the signal has stabilized or after pressing **Read**. The relevant buffer value is displayed and stored; the measurement icon disappears from the display.
 - To perform a 2-point calibration rinse the electrode with deionized water, place the electrode in the next calibration buffer and repeat the step above.
- 1 To finish the calibration and return to the sample measurement, press **Read**.
⇒ The offset value and the slope are then shown on the display for 3 seconds.
 - 2 To reject the calibration, press **Exit** before the meter returns automatically after 3 seconds to the measurement screen.

Note

With the 1-point calibration only the offset is adjusted. If the sensor was previously calibrated with multi-point calibration the previously stored slope will remain. Otherwise theoretical slope (-59.16 mV / pH) will be used.

5.1.4 Performing a 3-Point Calibration

- 1 Step 1: Perform the first and second calibration steps as described above in "Performing a 1 or 2-Point Calibration (page 9)".

- 2 Step 2: Rinse the electrode with deionized water.
- 3 Step 3: Place the electrode in the next calibration buffer and press **Cal**.
 - ⇒ The measurement icon appears.
The meter endpoints according to the preselected endpoint mode after the signal has stabilized or after pressing **Read**. The relevant buffer value is displayed and stored; the measurement icon disappears from the display.
The offset value and the slope are then shown on the display for 3 seconds.
- 4 To reject the calibration, press **Exit** before the meter returns automatically after 3 seconds to the measurement screen.

Note

- The use of a temperature sensor or an electrode with a built-in temperature sensor is recommended. If you use the **MTC** mode, you should enter the correct temperature value and keep all buffer and sample solutions at the set temperature.
- To ensure the most accurate pH readings, you should perform a calibration once a day.

5.2 Sample Measurement

5.2.1 Performing a pH Measurement

- To switch from mV mode to pH, press **Mode**.
- 1 Place the electrode in the sample and press **Read** to start the measurement:
⇒ The measurement icon appears on the display and the decimal point blinks. The display shows the pH value of the sample.
- 2 The automatic endpoint **A** is the default setting of the meter. When the signal has stabilized, the display freezes automatically, **/A** appears and the measurement icon disappears.
- or -
To manually endpoint a measurement, press **Read**. The display freezes and **/** appears.

Note

- By pressing and holding **Read**, you can switch between the auto and manual endpoint modes.
- Stability criterion for pH and mV measurement – The signal of the sensor input may not change by more than 0.15 mV in 6 seconds.

5.2.2 Performing a mV / Redox Measurement

- To switch from pH mode to mV, press **Mode**.
- For Redox measurements, connect a Redox sensor first.
- To perform a mV measurement, follow the same procedure as for pH measurement, see "Performing a pH Measurement (page 10)".

5.3 Temperature Measurement

For better accuracy, we recommend the use of either a built-in or a separate temperature probe. If a temperature probe is used, and the sample temperature are displayed.

Note

The meter accepts an NTC 30 kΩ temperature sensor.

5.3.1 Manual Temperature Compensation

If the meter does not detect a temperature probe, it automatically switches to the manual temperature compensation mode and **MTC** appears.

- 1 To set the **MTC** temperature and the temperature unit, press and hold **Mode/Setup** until the setup icon appears on the display and the **MTC** temperature blinks.
- 2 Use **▲** or **▼** to increase or decrease the value to the temperature of your sample. Press **Read** to confirm the setting.
⇒ The meter proceeds automatically to the temperature unit setting.
- 3 Press **Exit** to leave the setup without changing the temperature unit.
– or –
Use **▲** or **▼** to select between °C (Celsius) or °F (Fahrenheit). Press **Read** to confirm the setting.
⇒ The meter proceeds automatically to the calibration setting.
- 4 Press **Exit** to leave to setup. The default setting is 25 °C.

5.4 Using the Memory

5.4.1 Storing a Reading

FiveEasyPlus™ meter can store up to 99 endpointed results.

- Press **STO** when the measurement has endpointed.
⇒ **M01** indicates that one result has been stored.

Note

If you press **STO** when **M99** is displayed, **FUL** indicates the memory is full. To store further data you will have to clear the memory (see "Clearing the Memory (page 11)").

5.4.2 Recalling from Memory

- 1 Press and hold **RCL** to recall the stored values from memory when the current measurement has endpointed.
- 2 Press **▲** or **▼** to scroll through the stored results. **R01** to **R99** indicates which result is being displayed.
- 3 Press **Read** to exit.

5.4.3 Clearing the Memory

- 1 Continue pressing **▲** or **▼** to scroll through the stored results until **MRCL** appears.
- 2 Then press **Read**, **Clr** blinks.
- 3 Press **Read** again to confirm the deletion or press **Exit** to return to the measurement mode without deleting the data.

5.5 Printing

5.5.1 Connection and Configuration

A printer can be connected to the RS-232 interface of the instrument. The use of the **RS-P25**, **RS-P26**, **RS-P28** printer is recommended, because they recognize and adjusted to the correct parameters automatically.

In case a different printer is used, the following parameters need to be set:

Baud rate: 1200 bps
 Data bit: 8 bit
 Parity bit: None
 Stop bit: 1 bit

5.5.2 Printing a Measurement Results

If a printer is connected to the FEP20, a print-out is automatically generated after each endpointed measurement or calibration.

The format for the print-out following a pH measurement is:

EP, Value, Unit, Temp.
 AE 4.01pH 25.0°C MTC

The details for the second line are:

AE	4.01pH	25.0°C	MTC
			----- ATC or MTC
		-----	Temperature value
		-----	pH value
-----			AE for Auto-End-Point, ME for Manual-End-Point

The format for the print-out following a mV measurement is:

EP, Value, Unit, Temp.
 AE 930mV 25.0°C MTC

The print-out in case of an error message is:

EP, Value, Unit, Temp.
 Error2

5.5.3 Printing Calibration Results

The print-out for a 2-point calibration is:

Buffer1	4.01pH
mV1	178mV
Temp.1	25.0°C
Buffer2	7.00pH
mV2	-0mV
Temp.2	25.0°C
Slope	101%
Offset	-0mV
ATC/MTC	MTC

The print-out for a 3-point calibration is:

Buffer1	4.01pH
mV1	178mV
Temp.1	25.0°C
Buffer2	7.00pH
mV2	-0mV
Temp.2	25.0°C
Buffer3	9.21pH
mV3	-130mV
Temp.3	25.0°C
Slope	100%
Offset	-0mV
ATC/MTC	ATC

5.5.4 Printing from Memory

When scrolling through the memory (see "Recalling from Memory (page 11)") you can print the entry that is currently viewed by pressing and holding the **STO/RCL** key for 2 s.

The printout format is following:

Recall 01:
EP, Value, Unit, Temp.
ME 4.01pH 35.6°C ATC

5.6 Self-Diagnosis

- 1 Press and hold **Read** and **Cal** simultaneously until the meter displays the full screen.
 - ⇒ Each icon blinks one after the other. This way you may check whether all icons are correctly shown. The next step is to check that the keys are functioning correctly. This requires user interaction.
- 2 When **b** blinks, five icons are displayed. Press the five keys in any order.
 - ⇒ Each time you press a key an icon disappears from the screen.
- 3 Continue to press the other keys until all the icons have disappeared.
 - ⇒ When the self-diagnosis has been completed successfully, **PAS** appears. If self-diagnosis fails, error message **Err 1** appears (see "Error Messages (page 13)").

Note

You have to finish pressing all five keys within two minutes, otherwise **Err 1** appears and you will have to repeat the procedure.

5.7 Error Messages

Error 0	Memory access error	Reset to factory settings
Error 1	Self-diagnosis failed	Repeat the self-diagnosis procedure and make sure that you finish pressing all five keys within two minutes.
Error 2	Measured values out of range	Make sure that the electrode wetting cap has been removed and the electrode is properly connected and placed in the sample solution. If no electrode is connected, put the shorting plug in the socket.
Error 3	Measured buffer temperature out of range (5 ... 40 °C)	Keep the buffer temperature within the range for calibration
Error 4	Offset out of range	Make sure you have the correct buffer and that it is fresh. Clean or replace the electrode.
Error 5	Slope out of range	Make sure you have the correct buffer and that it is fresh. Clean or replace the electrode.
Error 6	Meter cannot recognize the buffer	Make sure you have the correct buffer and that it is fresh. Check that the buffer has not been used more than once during the calibration.

Error 9	The current data set has already been stored once	A measurement can only be stored once. Perform a new measurement to store a new data set.
---------	---------------------------------------------------	-------------------------------------------------------------------------------------------

6 Maintenance

6.1 Meter Maintenance

Never unscrew the two halves of the housing!

The Five series instruments do not require any maintenance other than occasional wipe with a damp cloth.

The housing is made of acrylonitrile butadiene styrene/polycarbonate (ABS/PC). This material is attacked by some organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK). Any spillage should be immediately wiped off.

6.2 Electrode Maintenance

Make sure the electrode is filled with electrolyte solution if applicable.

For maximum accuracy, any filling solution that may have "crept" and encrusted the outside of the electrode should be removed with deionized water.

Always store the electrode according to the manufacturer's instructions and do not allow it to dry out.

If the electrode slope falls rapidly, or if the response becomes sluggish, the following procedures may help. Try one of the following, depending the samples usually measured.

1. For fat or oil build-up, degrease the membrane with cotton wool soaked in either acetone or a soap solution.
2. If the sensor membrane has dried out, soak the tip of the electrode in 0.1 M HCl overnight.
3. If a protein build-up has occurred in the diaphragm, remove deposits by soaking the electrode in an HCl /pepsin solution (51350100).
4. If a silver sulfide contamination has occurred, remove deposits by soaking the electrode in a solution of thiourea (51350102).

Note

- After treatment a new calibration should be performed.
- Cleaning and filling solutions should be handled with the same care as that given to toxic or corrosive substances.

6.3 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

7 Accessories

Item	Order No.
Power supply	51302950
Electrode arm	51302951
Rubber covers (for electrode arm hole)	51302952
Shorting plug BNC	51302859
pH 4.01 buffer sachets, 30 x 20mL	51302069
pH 4.01 buffer solution, 6 x 250mL	51350018
pH 7.00 buffer sachets, 30 x 20mL	51302047
pH 7.00 buffer solution, 6 x 250mL	51350020
pH 9.21 buffer sachets, 30 x 20mL	51302070
pH 9.21 buffer solution, 6 x 250mL	51350022
pH 10.01 buffer sachets, 30 x 20mL	51302079
pH 10.01 buffer solution, 6 x 250mL	51350024
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/9.21)	51302068
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/10.01)	51302080
LE407, combination pH sensor, plastic, gel electrolyte	51340330
LE408, combination pH sensor, plastic, liquid electrolyte, refillable	51340347
LE409, combination pH sensor, glass, refillable	51340331
LE420, special pH sensor with PTFE junction, refillable	51340332
LE427, spear tip design, low maintenance	51340333
LE438, 3-in-1 pH sensor, plastic, ATC	51340242
LE501, Redox-sensor, glass	51340338
LE510, Redox-sensor, plastic	51340339
ATC probe, temperature sensor	51300164
HCl/Pepsin solution (removes protein contamination), 1 x 250 mL	51350100
Thiourea solution (remove silver sulfide contamination), 1 x 250 mL	51350102
Reactivation solution for pH electrodes, 1 x 25 mL	51350104

8 Specifications

FiveEasyPlus™ pH meter FEP20		
Measurement range	pH mode	pH 0.00...14.00
	mV mode	-1999...1999 mV
	Temperature	0 °C to 100 °C
Resolution	pH mode	0.01 pH
	mV mode	1 mV
	Temperature	0.1 °C
Limits of error	pH mode	± 0.01 pH
	mV mode	± 1 mV
	Temperature	± 0.5 °C
pH calibration	Up to 3 points	
Isopotential point	pH 7.00	
Calibration buffer	4 predefined groups	
Number of memory	99	
Power requirements	9 V DC, 0.4 W ± 10% For use with CSA certified (or equivalent approved) power source, which must have a limited circuit output.	
Size/weight	200 x 175 x 52 mm / 0.6 kg	
Display	Liquid crystal	
pH Input	BNC, impedance > 10e+12 Ω	
T Input	Cinch, NTC 30 kΩ	
Reference Input	2 mm banana socket	
IP rating	For indoor use only	
Ambient conditions	Temperature	5...40 °C
	Relative humidity	5 %...80 % (non-condensing)
	Height	up to 2000 m
	Installation category	II
	Pollution degree	2
Materials	Housing	ABS/PC reinforced
	Window	Membrane
	Keypad	Membrane

9 Appendix

9.1 Buffer Tables

FiveEasyPlus™ pH meter automatically corrects for the temperature dependence of the buffer pH using the values given in the tables.

Buffer group 1 (ref. 25 °C) METTLER TOLEDO US

5	7.09	4.00	10.25	1.67
10	7.06	4.00	10.18	1.67
15	7.04	4.00	10.12	1.67
20	7.02	4.00	10.06	1.68
25	7.00	4.00	10.01	1.68
30	6.99	4.01	9.97	1.68
35	6.98	4.02	9.93	1.69
40	6.97	4.03	9.89	1.69
45	6.97	4.04	9.86	1.70
50	6.97	4.06	9.83	1.71

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pH, conductivity & DO meters

EL20, EL2, FG4 (FEPx0, FEx0, ELx0, FGx, ELx -Series)

For additional types, see page type code

For optional equipment, see page accessories

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2006/95/EC Low voltage (LVD)

2004/108/EC Electromagnetic compatibility (EMC)

and that the standards have been applied.
und die Normen zur Anwendung gelangten.

Last two digits of the year in which the CE marking was affixed: 07

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2. T-H.

Rolf Truttmann
Head SBU pH Lab

References of standards for this declaration of conformity, or parts thereof:
Harmonized standards of Europe and Switzerland:

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Tested types

Portable equipment
FGx
ELx

Portable equipment
Series name for FiveGo
Series name for Education Line (Design is different)

Bench equipment
FEx0
ELx0

Bench equipment
Series name for FiveEasy
Series name for Education Line (Design is different)

Bench equipment
FEPx0

Bench equipment, software and with RS-232 port
Series name for FiveEasy Plus

x

Type of meter: 2 = ph, 3 = conductivity, 4 = DO
(Sensor and software are different)

Where x in the model designation may be any
number 0 to 9, denoting SELV/ELV secondary
circuits or minor mechanical differences.

Remarks

Bemerkungen:

Bench equipment: (0.6Kg / 200 x 175 x 52mm)
Rating: 9Vdc, 0.4W

Portable equipment: (0.18Kg / 169 x 82 x 36mm)
Rating: 6Vdc, 5mA powered with batteries:
4 x AAA / LR3 1.5V or
4 x NiMH 1.2V rechargeable

EC - DECLARATION OF CONFORMITY

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KD-Nr.: 51302969 D

Doku-Nr.: 20070006

Accessories

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Bench equipment:

FiveEasy™ Power supply, FiveEasy™ electrode arm
Rubber covers (for electrode arm hole), Shorting plug BNC

Portable equipment:

FiveGo™ electrode clip, Rubber clip covers (for electrode clip slot)
Education Line Battery cover, Wrist strap
FiveGo™ carry bag, IP54 electrode caps
Sealing kit (O-rings for battery case)

pH 4.01 buffer sachets, 30 x 20mL, pH 4.01 buffer solution, 6 x 250mL
pH 7.00 buffer sachets, 30 x 20mL, pH 7.00 buffer solution, 6 x 250mL
pH 9.21 buffer sachets, 30 x 20mL, pH 9.21 buffer solution, 6 x 250mL
pH 10.01 buffer sachets, 30 x 20mL
pH 10.01 buffer solution, 6 x 250mL
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/9.21)
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/10.01)

LE407, combination pH sensor, plastic, gel electrolyte

LE408, combination pH sensor, plastic, liquid electrolyte, refillable

LE409, combination pH sensor, glass, refillable

LE420, special pH sensor with PTFE junction, refillable

LE42/, spear tip design, low maintenance

LE438, 3-in-1 pH sensor, plastic, ATC

LE501, Redox-sensor, glass

LE510, Redox-sensor, plastic

ATC probe, temperature sensor

84 µS/cm standard solution, 500mL

1413 µS/cm standard solution sachets, 30 x 20 mL

1413 µS/cm standard solution, 250 mL

12.88 mS/cm standard solution sachets, 30 x 20 mL

12.88 mS/cm standard solution, 250 mL


LE703, conductivity sensor

LE 740, conductivity senso

To protect your product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of this product for years to come.

Please request full details about our attractive terms of service.

www.mt.com/phlab

For more information

Mettler-Toledo AG, Analytical
CH-8603 Schwerzenbach, Switzerland
Tel. +41 (0)44 806 77 11
Fax +41 (0)44 806 73 50
www.mt.com

Subject to technical changes.
© Mettler-Toledo AG 05/2012
51710598E