English

Operating instructions

Printer P25, P26, P28

Deutsch

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Drucker P25, P26, P28

Français

Mode d'emploi

Imprimante P25, P26, P28

Español

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Impresora P25, P26, P28

Italiano

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Stampante P25, P26, P28

English

Deutsch

Français

Español

Italiano

METTLER TOLEDO

Contents

Overview	4	6 Maintenance	14
1 Introduction	5	6.1 Changing the Paper Roll	15
1.1 Capabilities of the P25, P26, P28 P rinter.....	5	6.2 Replacing the Ribbon.....	15
1.2 Overview of Functions.....	5	6.3 Cleaning	15
2 Startup Procedure	5	6.4 Disposal	15
2.1 Safety Measures.....	5	7 Troubleshooting	16
2.2 Power Supply	6	8 Recommended Settings and Printer Models	17
2.3 Installing the Printer	6	9 Consumables and Spare Parts	19
3 Keystroke Functions	8	10 RS232C Interface	19
4 Menu	11	11 Technical Data	20
5 Printer Test	14	12 Declaration of Conformity	21

Overview

Controls

RS-P26

Connections

- 1 On/Off Switch
- 2 Pilot lamp – status indicator
- 3 Control panel
- 4 Power connection
- 5 Interface connector

Functions

Key	P25	P26	P28
	<ul style="list-style-type: none"> Paper feed Menu option (+ On Switch) 	<ul style="list-style-type: none"> Paper feed Scroll up menu option 	<ul style="list-style-type: none"> Paper feed Scroll up menu option
	—	<ul style="list-style-type: none"> Zeroing / taring the balance Abort – exit 	<ul style="list-style-type: none"> Zeroing / taring the balance Abort – exit
	—	<ul style="list-style-type: none"> Printout of date and time Scroll down menu option 	<ul style="list-style-type: none"> Printout of date and time Scroll down menu option
	—	<ul style="list-style-type: none"> Printing of stable weight values Opening the menu – select next menu option Save settings 	<ul style="list-style-type: none"> Printing of stable weight values Opening the menu – select next menu option Save settings
	—	—	<ul style="list-style-type: none"> Totaling application
	—	—	<ul style="list-style-type: none"> Statistics application

1 Introduction

1.1 Capabilities of the P25, P26, P28 Printer

This Printer is a simple to operate dot matrix printer for METTLER TOLEDO instruments. The print-out on normal paper meets the requirements of modern quality assurance systems (GLP, GMP, ISO 9001, etc.). These norms define the principles and user requirements related to measuring and test instruments. The Printers comply with the requirements of measurement traceability by precisely documenting the origin of data, as well as calibration procedures. It meets the requirement of archiving measurement results by providing a thermally stable and light resistant printout.

1.2 Overview of Functions

	RS-P25	USB-P25	RS-P26	RS-P28
Auto baudrate detection ¹⁾	✓	*	✓	✓
Balance feature "on/off" ²⁾	✓	*	✓	✓
USB Interface	–	✓	–	–
RS232 Interface	✓	–	✓	✓
Zero / Tare button	–	–	✓	✓
Date / Time	–	–	✓	✓
Statistics / Totaling	–	–	–	✓

* USB has standard settings

¹⁾ Auto baudrate detection: The Printer detects the settings of the connected device every time it is switches on/off by the key. Switch this feature off if the connected device has this feature not implemented (see chapter 8).

²⁾ Balance feature "on/off": This feature disables the keys on the Printer (see chapter 8).

2 Startup Procedure

2.1 Safety Measures

The Printer has been tested for the connections and intended purposes documented in these Operating Instructions. However, it is your responsibility to perform the necessary tests related to the methods and purposes of its intended use. Therefore, you should observe the following safety measures.

- The Printer must be used indoors only, and never in hazardous areas.
- The electronics of the Printer are protected against the ingress of splashing water and dust. However, since the printer housing is not watertight, the printer should not be used in the presence of liquids.
- Comply with cleaning instructions (see chapter 6.3), protect power plug against wet conditions.
- Never open the printer housing. If you have any problems with your printer, please contact your responsible METTLER TOLEDO service.

2.2 Power Supply

Your printer is supplied with an AC adapter with a country-specific power cable. The power supply is suitable for all line voltages in the range: 100 - 240 VAC, 50 - 60 Hz (for exact specifications, see chapter 9).

Attention

- First, check the local line voltage is in the range 100 - 240 VAC, 50 - 60 Hz and whether the power plug fits your local power supply connection. **If this is not the case, on no account connect the printer or the AC adapter to the power supply**, but contact the responsible METTLER TOLEDO dealer.
- **Only plug the adapter into a socket which is grounded.**

Important

- Before operating, check all cables for damage.
- Guide the cables so that they cannot become damaged or interfere with the weighing process!
- Take care that the AC adapter cannot come into contact with liquids!
- The power plug must be always accessible.

2.3 Installing the Printer


The Printer detects the settings of the Instrument. Make sure that you proceed as follow:

Connecting the printer to the instrument


- Attach the interconnect cable to the printer interface of the instrument and to the printer.
- Switch the instrument on

Switching the printer on

- Connect the Printer to the mains power.
- Use the main switch <I/O> to power on.
- RS232 interface only: The pilot lamp will flash until an instrument is connected. If the lamp is still flashing check chapter 7.

Inserting paper

- Remove paper cover (pull on the back) (1)
- Insert the paper rotary axis through the core of the paper roll (2).
- Feed paper through the slot in the printing unit in the direction of the arrow (3).
- Press and hold until enough paper has been fed (4).
- Lead paper through the paper cover slot
- Replace the paper cover (1).


Notes

- If, after the first printouts the connection settings (Bauderate, Bit, Handshake) of the balance change, the printer must be restarted (switch off / on) to start the auto bauderate detection again.
- Auto bauderate can be switched off.






3 Keystroke Functions

Printer P25

	 Press briefly	 Press and hold down
Operation Level	<p>  one line feed </p> <p>   Printer test with setting printout </p> <p>  On </p>	<p>  continuous feed (until release) </p> <p>   open menu </p> <p>  On </p>
Menu Level	<p>  next menu item </p>	<p>  store </p>

	 Press briefly	 Press and hold down
Operation Level	<div data-bbox="176 268 266 379">  zero / tare the balance </div> <div data-bbox="176 400 266 507">  print Date and Time </div> <div data-bbox="176 528 266 635">  one line feed </div> <div data-bbox="176 667 266 778">  print </div>	<div data-bbox="680 528 770 635">  continuous feed (until release) </div> <div data-bbox="680 667 770 778">  open menu </div>
Menu Level	<div data-bbox="165 863 613 1066" style="border: 1px dashed black; padding: 10px; display: flex; justify-content: space-around; align-items: center;"> <div data-bbox="176 863 266 986">  exit </div> <div data-bbox="281 863 370 986">  down </div> <div data-bbox="385 863 474 986">  up </div> <div data-bbox="489 863 579 986">  next </div> </div> <div data-bbox="680 863 770 986" style="margin-top: 20px;">  store </div>	

Printer P28

	 Press and hold down
<div style="display: flex; align-items: center;"> <div style="writing-mode: vertical-rl; transform: rotate(180deg); background-color: black; color: white; padding: 5px; font-weight: bold;">Operation Level</div> <div style="margin-left: 20px;"> <div style="margin-bottom: 10px;">  <p>zero / tare the balance</p> </div> <div style="margin-bottom: 10px;">  <p>print Date and Time</p> </div> <div style="margin-bottom: 10px;">  <p>one line feed</p> </div> <div>  <p>print</p> </div> </div> </div>	<div style="margin-bottom: 20px;">  <p>continuous feed (until release)</p> </div> <div>  <p>open menu</p> </div>
<div style="display: flex; align-items: center;"> <div style="writing-mode: vertical-rl; transform: rotate(180deg); background-color: black; color: white; padding: 5px; font-weight: bold;">Menu Level</div> <div style="margin-left: 20px;"> <div style="display: flex; justify-content: space-around; align-items: center; border: 1px dashed gray; padding: 10px;"> <div style="text-align: center;">  <p>exit</p> </div> <div style="text-align: center;">  <p>down</p> </div> <div style="text-align: center;">  <p>up</p> </div> <div style="text-align: center;">  <p>next</p> </div> </div> </div> </div>	<div style="margin-bottom: 20px;">  <p>store</p> </div>
<div style="display: flex; align-items: center;"> <div style="writing-mode: vertical-rl; transform: rotate(180deg); background-color: black; color: white; padding: 5px; font-weight: bold;">Function Level</div> <div style="margin-left: 20px;"> <div style="margin-bottom: 10px;">  <p>Start totaling application</p> </div> <div>  <p>Start statistics application</p> </div> </div> </div>	<div style="margin-bottom: 20px;">  <p>Conclude totaling application</p> </div> <div>  <p>Conclude statistics application</p> </div>

4 Menu


Printer P25

 Press briefly	 Press and hold down	Printout
<div style="border: 1px solid black; padding: 5px;"> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Settings</div> <div style="padding: 2px; text-align: center;">Reset</div> <div style="text-align: center;">↓</div> <div style="border: 1px dashed gray; padding: 2px; text-align: center;"> RS232 interface only </div> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Auto Bauderate</div> <div style="padding: 2px; text-align: center;">On</div> <div style="text-align: center;">↓</div> <div style="padding: 2px; text-align: center;">Auto Bauderate Off</div> <div style="text-align: center;">↓</div> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Char Set</div> <div style="padding: 2px; text-align: center;">IBM / DOS</div> <div style="text-align: center;">↓</div> <div style="padding: 2px; text-align: center;">IBM / DOS Cyrillic</div> <div style="text-align: center;">↓</div> <div style="padding: 2px; text-align: center;">IBM / DOS Japanese</div> <div style="text-align: center;">↓</div> <div style="padding: 2px; text-align: center;">ANSI / WIN Latin 1</div> <div style="text-align: center;">↓</div> <div style="padding: 2px; text-align: center;">ANSI / WIN Cyrillic</div> <div style="text-align: center;">↓</div> <div style="padding: 2px; text-align: center;">ANSI / WIN Japanese</div> <div style="text-align: center;">↓</div> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Balance Feature</div> <div style="padding: 2px; text-align: center;">Enable</div> <div style="text-align: center;">↓</div> <div style="padding: 2px; text-align: center;">Balance Feature Disable</div> </div>	<p style="text-align: center;">Reset the printer</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p> <p style="text-align: center;">Store</p>	<p>Settings Reset</p> <p>----- Reset done ----- Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS</p> <p>Auto Bauderate: On -----STORED----- ...<i>"current settings"</i></p> <p>Auto Bauderate: Off -----STORED----- ...<i>"current settings"</i></p> <p>IBM/DOS -----STORED----- ...<i>"current settings"</i></p> <p>IBM/DOS Cyrillic -----STORED----- ...<i>"current settings"</i></p> <p>IBM/DOS Japanese -----STORED----- ...<i>"current settings"</i></p> <p>ANSI/WIN Latin 1 -----STORED----- ...<i>"current settings"</i></p> <p>ANSI/WIN Cyrillic -----STORED----- ...<i>"current settings"</i></p> <p>ANSI/WIN Japanese -----STORED----- ...<i>"current settings"</i></p> <p>Balance Feature: Enable -----STORED----- ...<i>"current settings"</i></p> <p>Balance Feature: Disable -----STORED----- ...<i>"current settings"</i></p>

English

Menu Option
Factory Setting

Printer P26 / P28

Examples if balance display is used:

Settings	
Reset	rESEt
List	LlSt
Header	
*No Header	no Hdr
Date/Time	d t Hdr
Date/Time/Balance	d t b Hdr
Date Day**	
18	dRy 18
Date Month**	
07	Month 07
Date Year**	
2008	YEAR-2008
Time Minute**	
28	Min. 28
Time Hour**	
14	Hour 14
Date Format***	
DD.MMM.YYYY	ddMMY
MMM DD YYYY	MM dd Y
DD.MM.YYYY	ddMMY
MM/DD/YYYY	MM/ddYY
YYYY-MM-DD	Y-MM-dd
Time Format***	
24:MM:SS	24:MM:SS
12:MM:SS	12:MM:SS
24.MM.SS	24MM.SS
12.MM.SS	12MM.SS
Baud rate	
*9600	bd9600
19200	bd.19200
38400	bd38400
600	bd600
1200	bd.1200
2400	bd2400
4800	bd4800
Bit/Parity	
*8 Bit No	8b-no
7 Bit Odd	7b-odd
7 Bit Even	7b-E
7 Bit No	7b-no
Handshake	
*Xon/Xoff	HS oFF
HS Soft	HS 5oFt
HS Hard	HS HRRd
Char Set	
*IBM/DOS	d05
IBM/DOS Cyrillic	d05 C4r
IBM/DOS Japanese	d05 JRP
ANSI/Win Latin 1	11 in LRt
ANSI/Win Cyrillic	11 in C4r
ANSI/Win Japanese	11 in JRP
Tare/Zero	
*Zero	ZEr0
Tare	tAR-E



Menu Option
Factory Setting

* factory settings



** instruments with built-in clock: not applicable, built-in clock has priority

*** instruments with built-in clock: format of printer has priority

Example for setting the day (P26 / P28 only)

	Press briefly 	Press and hold down 
P26 / P28	Balance Display	
	 * MENU ▼ * rESEt	Open the menu: → Hold down the «  » key until "reset" appears.
	 * dAY 18	→ Select the menu option "Day" with the «  » key.
	 * dAY 21	→ Select the "Day" – scroll up with the «  » key.
		– scroll down with the «  » key.
	 * StorEd	Save the settings: → Hold down the «  » key until "Stored" appears.

Applications (P28 only)

 Press briefly	 Press and hold down
Statistics Application ¹⁾	Totaling Application ²⁾
	
→ Start "Statistics application" (LED on).	→ Start "Totaling application" (LED on).
	
→ Add a new sample and press «M+».	→ Add a new sample and press «M+».
	
→ If needed print subtotal press the «  » key.	→ If needed print subtotal press the «  » key.
	
→ Print and end "Statistics application" (LED off).	→ Print and end "Totaling application" (LED off).

Important notes

- Unit change during the application is not permitted.
 - Maximal 999 sample.
- ¹⁾ • Sample weight must be in the range of 70% - 130% of the current average value.
- Minimum dynamic weight changed must be 100 d.
- ²⁾ Minimum weights + 1 d.

5 Printer Test

The Printer is equipped with a self-test that automatically prints out the character set of the printer. (standard IBM/DOS).

Starting the test

- Press the key «↑» and power on the Printer.
The printer prints out the current settings.

Terminating the test

- Power off the Printer.

6 Maintenance


Under normal conditions the Printer requires practically no maintenance. The METTLER TOLEDO service department is available if servicing is required. Please ask your responsible METTLER TOLEDO dealer for details.

6.1 Changing the Paper Roll

Insertion of a new roll of paper is described in chapter 2.3. See chapter 9 (Consumables) for the paper roll's order number.

Note: You may obtain the paper rotary axis from the old paper roll.

6.2 Replacing the Ribbon

- Remove the paper cover (1).
- Pull paper out of printing unit (2).

- To remove the ribbon cartridge, press in direction of arrow (3).
- Insert the new ribbon, and if necessary, take up slack with tensioning wheel (4). See chapter 9 (Consumables) for the ribbon cartridge order number.
- Insertion of a new roll of paper is described in chapter 2.3.

6.3 Cleaning

Since the printer housing is made of top grade, resistant materials, all commercially available mild cleaning agents may be used (see safety precautions in chapter 2.1).

6.4 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

If this device is passed on to other parties (for private or professional use), compliance with this regulation will be required as well.

Thank you for your contribution to environmental protection.

7 Troubleshooting

Problem / Message	Possible cause	Correction
Printout unreadable	<ul style="list-style-type: none"> Ribbon worn or entangled. Service life of printing unit at an end. 	<ul style="list-style-type: none"> Change ribbon cartridge or tense the ribbon properly. Please contact METTLER TOLEDO service.
Pilot lamp does not light	<ul style="list-style-type: none"> No power. Power supply fault. 	<ul style="list-style-type: none"> Switch printer on. Please contact METTLER TOLEDO service.
Pilot lamp flashes ---Bauderate not found--- *Please connect device and switch the printer off/on. *Alternatively you can change the settings of the printer or the connected device. Printer settings: Baudrate 9600 Bit/Parity 8 Bit None Handshake Off	<ul style="list-style-type: none"> No connection printer device. Printer did not find the settings of the connected instrument. 	<ul style="list-style-type: none"> Make sure that the printer and the instrument is connected. Switch the printer off/on. Alternatively change the settings of the printer or the connected instrument. (see chapter 8)
-----Timeout -----	<ul style="list-style-type: none"> Data transfer has not completed in the last 45 s. 	<ul style="list-style-type: none"> Move balance to more stable location. Check the balance.
----- Over Load -----	<ul style="list-style-type: none"> Balance in overload range. 	<ul style="list-style-type: none"> Check the weight of sample.
	<ul style="list-style-type: none"> Zeroing out of the zero range. 	<ul style="list-style-type: none"> Change the zero range of the balance. Change the setting of the printer from zero to tare.
----- Under Load -----	<ul style="list-style-type: none"> Balance in underload range. Weighing pan is not in position. 	<ul style="list-style-type: none"> Check the balance.
-- Transmission Error --	<ul style="list-style-type: none"> Instrument does not understand printer instruction. 	<ul style="list-style-type: none"> Check interface parameter settings on printer and instrument.
---- Not Executable ----	<ul style="list-style-type: none"> Instrument can not execute printer instruction. 	<ul style="list-style-type: none"> Check the balance.
--- Date/Time Not Set --- ----- Memory Error -----	<ul style="list-style-type: none"> Printer error messages. 	<ul style="list-style-type: none"> Check the instrument. Reset the printer. Set date/time. Back-up battery depleted Please contact METTLER TOLEDO service.
----- Memory full -----	<ul style="list-style-type: none"> Statistics or Totaling Memory is full. 	<ul style="list-style-type: none"> Start a new Statistics or Totaling application.
-- Load/Unload Weight --	<ul style="list-style-type: none"> No sample or the same sample is on the weighing pan. 	<ul style="list-style-type: none"> Put a new sample or remove it and put it again on the pan.

Problem / Message	Possible cause	Correction
Unit change not allowed !	<ul style="list-style-type: none"> Changed unit during Statistics or Totaling application. 	<ul style="list-style-type: none"> Reset the balance. Change the unit back to the first one
-- Out of range --	<ul style="list-style-type: none"> Sample weight is out of range 	<ul style="list-style-type: none"> Sample weight must be in the range of 70%–130% of the current average value.


8 Recommended Settings and Printer Models

Connected Device	Device Settings	Printer Settings	Compatible Printer Models
MS, ML, JP, JS, PHS, PHL, AB-S/FACT, AB-S, AB-L, PB-S/FACT, PB-S, PB-L, AL, PL, PL-S, PL-L, JB-G/FACT, JB-C/FACT, JB-L-G, JB-L-C, JL-G, JL-C, JL-G/L, EL, AB-S/PH, PB-S/PH	Printer Mode	Factory settings	RS-P25, RS-P26, RS-P28
PG-S, PG-S/PH	Printer Mode	Factory settings	RS-P25
SB	Printer Mode	Factory settings	RS-P26
XP, XS, XA	Printer Mode	Factory settings	RS-P25
Katakana	Printer Mode and language ANSI/WIN Katakana	Factory settings	
Russian	Printer Mode and language ANSI/WIN Russian	CharSet IBM/DOS Cyrillic	
Chinese, Japanese	Printer Mode	Factory settings	
MJ33, HB43-S	Printer Mode	Factory settings	RS-P25
Nihongo	Printer Mode and language Nihongo	Factory settings	
Russian	Printer Mode and language Russian	CharSet IBM/DOS Cyrillic	

Connected Device	Device Settings	Printer Settings	Compatible Printer Models
HG63, HR83	Select on LocalCAN Switch-box: 1 Printer, Baudrate: 9600; Bit/Parity: 8 Bit No, Handshake: Xon/Xoff	*Reset Auto baudrate off, Balancefeature off	RS-P25
Russian	Select on LocalCAN Switch-box: 1 Printer, Baudrate: 9600; Bit/Parity: 8 Bit No, Handshake: Xon/Xoff and language Russian	*Reset Auto baudrate off, Balancefeature off, CharSet IBM/DOS Cyrilic	
QB5	Factory settings	Factory settings	RS-P25
SQC XP	Printer Mode	Factory settings	RS-P25
SQC16	Printer Mode	*Reset, Auto baudrate off	RS-P25
T50, T70, T90	USB compact printer	Factory settings	USB-P25
V20, V30, C20, C30	USB compact printer	Factory settings	USB-P25
DL15, DL22, DL28	RS-P42	*Reset Auto baudrate off, Balancefeature off	RS-P25 (Connect with RS9 F/F Cable Order number: 51190362)
DL31, DL38, DL32, DL39	GA42	Auto baudrate off, Balancefeature off, Baudrate: 1200; Bit/Parity: 8 Bit No, Handshake: None	RS-P26 (Connect with RS9 F/F Cable Order number: 51190362)
DE40, DE45, DE51, RE40D, RE50, DR40, DR45	GA42	Auto baudrate off, Balancefeature off, Baudrate: 1200; Bit/Parity: 8 Bit No, Handshake: Hardware	RS-P26 Connect with RS9 F/F Cable (Order number: 51190362)
Refracto 30PX/GS, Densito 30PX	PrinterMode and Baudrate 9600	Auto baudrate off, Balancefeature off, Baudrate: 9600; Bit/Parity: 8 Bit, Handshake: Hardware	RS-P26
SevenEasy S20, S30, V 1.30 or higher	Factory settings	*Reset Auto baudrate off, Balancefeature off	RS-P26, RS-P25
SevenMulti S40, S47, S50, S70, S80, V 1.1 or higher	Baudrate: 9600; Bit/Parity: 8 Bit No, Handshake: Xon/Xoff	*Reset, Auto baudrate off, Balancefeature off	RS-P25, RS-P26
SevenGo pro SG6, SG7, SG8	Factory settings of the RS232 Infrared device	Auto baudrate off Bau- drate: 1200; Bit/Parity: 8 Bit No, Handshake: off	RS-P26 Datatransfer of SevenGo pro to the printer occurs via the Infrared-RS232 adapter (Order number: 51302333)
DM40, DM45, DM50, RM40, RM50	Printer Mode	Factory settings	USB-P25
MP50, MP70, MP80	Printer Mode	Factory settings	USB-P25

* The Printer detects automatically the settings of the connected device. Therefore disconnect the Printer of the connected device and reset the Printer to set the factory settings again: Baudrate: 9600; Bit/Parity: 8 Bit No, Handshake: Xon/Xoff

9 Consumables and Spare Parts


Pos.	Item	Part number	Delivery quantity
1	Paper cover	12120734	1 pcs
2	Paper roll 57.5 x ø 50 mm	00072456	5 pcs
	Paper roll, self-adhesive	11600388	1 pcs
3	Paper rotary axis	12120735	1 pcs
4	Ribbon cartridge, black	00065975	2 pcs
-	RS9 M/F Cable, 1.5 m	11101051	1 pcs
-	RS9 F/F Cable, 1.5 m	51190362	1 pcs
-	USB Cable	12130716	1 pcs
5	AC/DC adapter (without power cable) 100 - 240 VAC, 50 - 60 Hz, 12 VDC 2.5 A	11107909	1 pcs
6	Power cable CH	00087920	1 pcs
	Power cable EU	00087925	1 pcs
	Power cable US	00088668	1 pcs
	Power cable IT	00087457	1 pcs
	Power cable DK	00087452	1 pcs
	Power cable GB	00089405	1 pcs
	Power cable AU	00088751	1 pcs
	Power cable ZA	00089728	1 pcs
	Power cable BR	30015268	1 pcs
	Power cable JL	00225297	1 pcs
	Power cable IN	11600569	1 pcs
	Power cable JP	11107881	1 pcs
	Power cable TH, PE	11107880	1 pcs

English

10 RS232C Interface

RS-P25, RS-P26 and RS-P28 printers are equipped with an RS232C interface to connect METTLER TOLEDO instruments.

- 9-pin male connector
- Matching to other device (transmission parameters) see chapter 4.

11 Technical Data

Power supply	AC/DC adapter: 12 VDC, 2.5 A Primary: 100 - 240 VAC, -50 - 60 Hz Secondary: 12 VDC, 2.5 A (with electronic overload protection)
Cable to AC/DC adapter	3-core, with country-specific plug
Power line voltage	12 VDC
Power consumption	Max. 1 A
Altitude	up to 4000 m
Ambient temperature	0 °C – 40 °C
Relative humidity	max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, noncondensing
Overvoltage category	Class II
Pollution degree	2
Printing unit	Dot matrix printer 5x7, 24-character line length
Printing speed	2.3 lines per second
Ribbon cartridge	Exchangeable, black
Paper roll	Standard paper 57.5 mm x ø 50 to 60 mm, integrated in housing, commercial size
Print quality	Light resistant and thermally stable printing (GLP, GMP, ISO 9001)
Interface:	RS-P25, RS-P26, RS-P28 USB-P25
	RS232C USB
Printer dimensions	(WxDxH) 120 mm x 203 mm x 73 mm
Packaging dimensions	(WxDxH) 255 mm x 205 mm x 150 mm
Net weight	740 g (paper roll incl.)
Compatibility limits	<ul style="list-style-type: none"> • Function Date / Time: not applicable on balances with built-in clock (built-in clock has priority) • Function Zeroing / Taring: not applicable on all models of balances

12 Declaration of Conformity

The undersigned declare on behalf of
Mettler-Toledo AG
Im Langacher
CH-8606 Greifensee

that the Printer **METTLER TOLEDO P25, P26** and **P28** to which this declaration relates (serial number specified on the product) is in compliance with the below mentioned EC Directives (including all amendments)

2006/95/EC Low Voltage Directive

2004/108/EC Electromagnetic compatibility

and that the following standards have been applied

IEC/EN61010-1:2001,

IEC/EN61326-1:2005 / EN61326-1:2006 (class B)


IEC/EN61326-1:2005 / EN61326-1:2006 (industrial environments)

for Canada, USA and Australia

CAN/CSA-C22.2 No.61010-1-04, UL Std. No.61010A-1 (2nd Edition), FCC, Part 15, class A, AS/NZS CISPR 11, AS/NZS 61000.4.3

Greifensee, 09.03.2012

Mettler-Toledo AG
Laboratory & Weighing Technologies

René Lenggenhager
General Manager

Marcel Strotz
Manager SBU LAB Basic Weighing

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
4 Menü

Drucker P25

 Kurz drücken	 Lange drücken	Druckausgabe
<div style="border: 1px solid black; padding: 5px;"> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Settings</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Reset</div> <div style="text-align: center;">↓</div> <div style="border: 1px dashed gray; padding: 5px; margin: 5px 0;"> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Auto Bauderate</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">On</div> <div style="text-align: center;">↓</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Auto Bauderate Off</div> </div> <div style="text-align: center;">↓</div> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Char Set</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS</div> <div style="text-align: center;">↓</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS Cyrillic</div> <div style="text-align: center;">↓</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS Japanese</div> <div style="text-align: center;">↓</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Latin 1</div> <div style="text-align: center;">↓</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Cyrillic</div> <div style="text-align: center;">↓</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Japanese</div> <div style="text-align: center;">↓</div> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Balance Feature</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Enable</div> <div style="text-align: center;">↓</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature Disable</div> </div>	<p style="text-align: center;">Drucker rückstellen (Reset)</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p> <p>speichern</p>	<pre> Settings Reset ---- Reset done ---- Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS Auto Bauderate: On -----STORED----- ... "current settings" Auto Bauderate: Off -----STORED----- ... "current settings" IBM/DOS -----STORED----- ... "current settings" IBM/DOS Cyrillic -----STORED----- ... "current settings" IBM/DOS Japanese -----STORED----- ... "current settings" ANSI/WIN Latin 1 -----STORED----- ... "current settings" ANSI/WIN Cyrillic -----STORED----- ... "current settings" ANSI/WIN Japanese -----STORED----- ... "current settings" Balance Feature: Enable -----STORED----- ... "current settings" Balance Feature: Disable -----STORED----- ... "current settings" </pre>
		<div style="background-color: black; color: white; padding: 2px;">Menüpunkt</div> <div style="border: 1px solid black; padding: 2px;">Werkseinstellung</div>

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Bit/Parität

*8 Bit No *8b-no*
7 Bit Odd *7b-odd*
7 Bit Even *7b-E*
7 Bit No *7b-no*

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
4 Menu

Imprimante P25

 Appuyez brièvement	 Maintenez appuyé	Message imprimé
<div style="border: 1px solid black; padding: 5px;"> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Settings</div> <div style="padding: 2px; text-align: center;">Reset</div> </div>	<p style="text-align: center;">réinitialiser l'imprimante</p>	<pre>Settings Reset ---- Reset done ---- Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Auto Bauderate</div> <div style="padding: 2px; text-align: center;">On</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>Auto Bauderate: On -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Auto Bauderate</div> <div style="padding: 2px; text-align: center;">Off</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>Auto Bauderate: Off -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Char Set</div> <div style="padding: 2px; text-align: center;">IBM / DOS</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>IBM/DOS -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="padding: 2px; text-align: center;">IBM / DOS Cyrillic</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>IBM/DOS Cyrillic -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="padding: 2px; text-align: center;">IBM / DOS Japanese</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>IBM/DOS Japanese -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="padding: 2px; text-align: center;">ANSI / WIN Latin 1</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>ANSI/WIN Latin 1 -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="padding: 2px; text-align: center;">ANSI / WIN Cyrillic</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>ANSI/WIN Cyrillic -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="padding: 2px; text-align: center;">ANSI / WIN Japanese</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>ANSI/WIN Japanese -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="background-color: black; color: white; padding: 2px; text-align: center;">Balance Feature</div> <div style="padding: 2px; text-align: center;">Enable</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>Balance Feature: Enable -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 5px;"> <div style="padding: 2px; text-align: center;">Balance Feature Disable</div> </div>	<p style="text-align: center;">enregistrement</p>	<pre>Balance Feature: Disable -----STORED----- ... "current settings"</pre>

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4 Menu

Stampante P25

 Premere brevemente	 Premere e mantenere premuto	Stampa
<div style="border: 1px solid black; padding: 5px;"> <p style="text-align: center;">Settings</p> <p style="text-align: center;">Reset</p> <hr/> <div style="border: 1px dashed gray; padding: 2px; text-align: center; font-size: small;">Solo per l'interfaccia RS232</div> <p style="text-align: center;">Auto Bauderate</p> <p style="text-align: center;">On</p> <hr/> <p style="text-align: center;">Auto Bauderate</p> <p style="text-align: center;">Off</p> <hr/> <p style="text-align: center;">Char Set</p> <p style="text-align: center;">IBM / DOS</p> <hr/> <p style="text-align: center;">IBM / DOS</p> <p style="text-align: center;">Cyrillic</p> <hr/> <p style="text-align: center;">IBM / DOS</p> <p style="text-align: center;">Japanese</p> <hr/> <p style="text-align: center;">ANSI / WIN</p> <p style="text-align: center;">Latin 1</p> <hr/> <p style="text-align: center;">ANSI / WIN</p> <p style="text-align: center;">Cyrillic</p> <hr/> <p style="text-align: center;">ANSI / WIN</p> <p style="text-align: center;">Japanese</p> <hr/> <p style="text-align: center;">Balance Feature</p> <p style="text-align: center;">Enable</p> <hr/> <p style="text-align: center;">Balance Feature</p> <p style="text-align: center;">Disable</p> </div>	<p>Resettare la stampante</p> <p>memorizza</p> <p>memorizza</p> <p>memorizza</p> <p>memorizza</p> <p>memorizza</p> <p>memorizza</p> <p>memorizza</p> <p>memorizza</p> <p>memorizza</p> <p>memorizza</p>	<pre> Settings Reset ----- Reset done ----- Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS Auto Bauderate: On -----STORED----- ... "current settings" Auto Bauderate: Off -----STORED----- ... "current settings" IBM/DOS -----STORED----- ... "current settings" IBM/DOS Cyrillic -----STORED----- ... "current settings" IBM/DOS Japanese -----STORED----- ... "current settings" ANSI/WIN Latin 1 -----STORED----- ... "current settings" ANSI/WIN Cyrillic -----STORED----- ... "current settings" ANSI/WIN Japanese -----STORED----- ... "current settings" Balance Feature: Enable -----STORED----- ... "current settings" Balance Feature: Disable -----STORED----- ... "current settings" </pre>

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© Mettler-Toledo AG 03/2012
Mettler-Toledo AG Laboratory & Weighing Technologies