
Operating Instructions
Bedienungsanleitung
Mode d'emploi
Instrucciones de manejo

METTLER **TOLEDO**

PE360 DeltaRange
PE3600 DeltaRange

NOTE (for users in the U.S.A.)

This equipment generates and uses radio frequency energy and if not installed and used properly, that is, in strict accordance with the manufacturer's instructions, may cause interference to radio and television reception. It has been type tested and found to comply with the limits for a Class B computing device in accordance with the specifications in Subpart J of Part 15 of FCC Rules, which are designed to provide reasonable protection against such interference in a residential installation. However, there is no guarantee that interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient the receiving antenna.
- Relocate the equipment with respect to the receiver.
- Move the equipment away from the receiver.
- Plug the equipment into a different outlet so that equipment and receiver are on different branch circuits.

If necessary, the user should consult the dealer or an experienced radio/television technician for additional suggestions. The user may find the following booklet prepared by the Federal Communications Commission helpful:

"How to Identify and Resolve Radio-TV Interference Problems".

This booklet is available from the U.S. Government Printing Office, Washington, DC 20402, Stock No. 004-000-00345-4.

OPERATING INSTRUCTIONS

Page 2

English

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Seite 20

Deutsch

MODE D'EMPLOI

Page 38

Français

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Página 56

Español

SUMMARY OF MODELS



These Operating Instructions apply for the balances in the PE model line presented here. Operation of these balances is identical. The differences are in the weighing range, readability and size of the weighing pan.

PE360 DeltaRange

Weighing range	360 g
– DeltaRange	60 g
Readability	0.01 g
– in the DeltaRange	0.001 g

PE3600 DeltaRange


Weighing range	3600 g
– DeltaRange	600 g
Readability	0.1 g
– in the DeltaRange	0.01 g

The two balances, the PE360 DeltaRange and the PE3600 DeltaRange, allow you to weigh very easily in two different ranges: either the coarse range, or the fine range (DeltaRange) whose readability is 10 times greater and which can be shifted to any place within the entire balance weighing range.

These PE models are equipped with:

- **Nonmetric unit-of-weight displays**, which make it possible to indicate weights not only in metric units of weight but also in several nonmetric units. For the sake of simplicity, all operating procedures, information and technical data in these Operating Instructions are given in grams. However, the same also applies to nonmetric units of weight which can be selected by the user.
- **Autocalibration**, a semiautomatic calibration device used to adjust the weighing accuracy.
- **METTLER TOLEDO DeltaDisplay**, a special display mode which is useful when weighing-in.
- **Stability detector**, a signal indicating unstable weighing results.
- **OFF display**, a signal indicating that a power failure has occurred earlier.
- **Overload indication**, which tells the operator that the balance has exceeded its weighing range.
- **Underload indication**, which tells the operator that the balance has not reached its minimum load.
- **Plastic draft shield**, which prevents unstable weight indications caused by drafts (supplied as standards equipment in PE360 only).

Operating elements and connections

- (1) Single control bar
- (2) Connection socket of data output
- (3) Connection socket for
GE305 Application Input Device,
GE310 Terminal or GE48 Auxiliary display
- (4) Fuse holder
(with spare fuse)
- (5) Power-line connection socket

Items (2) and (3) are only installed in the balance if the instrument is equipped with a Data Output (available on request).

TABLE OF CONTENTS

Page

Preliminary steps

- Checking the voltage setting 4
- Installing the draft shield, setting up the balance 5

Operation

- Switching the balance on or off 6
- Calibrating the balance 7
- Selecting the weight unit 8
- Determining weight or taring 9
- Weighing-in using DeltaDisplay 10
- Making subtractive weighings, weight comparisons 11

Further capabilities

- Level indicator Field Installation Kit;
glass draft shield, GD hanger 12
- Application Packages; Data Output 13

Peripheral instruments

14

Accessories

15

Maintenance

- Installing the in-use cover and changing the fuse 16
- Cleaning and calibration 17

What's wrong if ...

18

Technical data

19

PRELIMINARY STEPS: Checking and correcting the voltage setting


The voltage set in the balance must agree with the power-line voltage.

Before putting the balance into operation for the first time, the voltage setting should be checked and, if needed, corrected.

The voltage set at the factory is listed on a label on the power connection socket (5).

If this voltage setting corresponds to your local power-line voltage, go on the next page.

If this voltage setting does not correspond to your local power-line voltage, or if the label is missing, the balance housing must be opened, the voltage set on the voltage selector must be checked, and if necessary, changed.


Opening the housing

Before doing so, make sure that the power-line cable is disconnected!

When the power-line cable is connected, the balance electronics are powered, even though the display does not light up!

- If the weighing pan (6) and pan support (7) are installed, remove them.
- Remove screw (8).
- Carefully remove housing cover (10) with in-use cover, if installed.


Voltage selector next to weighing cell

- Check whether the voltage selector jumper (11) is plugged onto the pin that is labeled with your local power-line voltage.
- If necessary, change position of voltage selector and plug it onto the correct pin.

Admissible power-line voltages:

95 V,	105 V	{ Install 125 mA slow
110 V,	120 V	{ blowing microfuse
190 V,	210 V	{ Install 63 mA slow
220 V,	240 V	{ blowing microfuse

Note: When switching the operating voltage from 95...120 V to 190...240 V (or vice versa), the microfuse must be changed. The two jumper plugs (12) must not be changed.

Closing the balance housing

- Carefully set the housing cover (10) down onto the balance.
- Insert screw (8) and tighten.

Installing the weighing pan

- If still in place, remove protective coating from weighing pan.
- Place the pan support (7) onto the conical peg (13).
- Place the weighing pan (6) on the pan support.


Connecting the power cable (14): see next page under "Balance location".

PRELIMINARY STEPS: Installing the draft shield and setting up the balance

Installing the pan support and weighing pan is always done in the same manner, regardless of whether a draft shield and/or in-use cover is used or not. With the PE360, the draft shield provided should be installed, so as to achieve measurement values that are not affected by drafts.

Every PE balance provides reliable measuring values, even under less than favorable ambient conditions.


However, it is preferable to set the balance up at a stable location, with as little draft as possible and that is not in direct sunlight.

Installing the draft shield (PE360)

The base (19) is installed in the factory. If it was subsequently removed, proceed as follows:


- Remove weighing pan and pan support.
- Place base (19) on the balance and turn it clockwise until it stops.
- Place pan support and weighing pan on balance.

Level indicator field installation kit

- Depending on the height of the draft shield required, place low or high middle section (20) on the base.
- Place cover (21) with the two snaps on the middle section.

If wished, a glass draft shield can be ordered for the PE360, see ACCESSORIES.


Balance location

- Stable base, that is vibration-free and horizontal.
- No large temperature fluctuations.
- Avoid direct sunlight.
- Find a location free of drafts.
- At the selected location, connect power cable (14).

OPERATION: How to switch the balance on or off


The accuracy of an electronic measuring instrument is more constant if the operating temperature remains the same. For this reason, these PE balances are built so that when the instrument is switched off using the single control bar (1), only the display is switched off. The balance electronic components remain switched on as long as the power cable is connected (stand-by operation).

It is also important that the balance display always be switched on (by pressing the control bar) with the pan empty.

Switching the balance on


- With the pan empty, briefly press the single control bar (1); all display segments light up for several seconds:

This allows you to make a check of the display.


Then the display indicates zero, with as many decimal places as the balance model readability permits, e.g.:

0.000 g with the PE360 or
0.00 g with the PE3600.

Switching off the display

- Lift single control bar (1); the display is switched off.

OFF Display


If there is a power interruption at any time, "OFF" will appear in the display as soon as power is restored. In this case, the single control bar (1) must be pressed briefly.

OPERATION: How to calibrate the balance

Your balance must be calibrated every time it is moved to a new location, i.e., the correct weight display must be adjusted.

Only then will your weighings be accurate. The balance must be left connected to the power source for 30 minutes before calibrating it (warm-up time). The plastic draft shield should be on the PE360.

Depending on the balance model, certain test weights are needed for calibration. These weights are listed under "ACCESSORIES, optional" for each balance model. The test weights should always be handled carefully and, as much as possible, should not be touched with your bare hands. If your balance is equipped with a leveling kit, it should be leveled before it is calibrated.


Checking calibration

- Press the single control bar (1); zero appears in the display.
- Place test weight on pan and read the display after the ring symbol (15) (stability detector) has been blanked out.

If the balance displays the weight off the test weight accurately to the last decimal place, calibration is okay. If not, the balance must be calibrated:


Correct calibration as described in the middle paragraph.

Correcting the calibration

- Press the control bar (1) until the display indicates "-----", then release the control bar.
- "- CAL -" appears in the display.
- Place test weight on pan, i.e.,:
 - 100 g for the PE360,
 - 1000 g for the PE3600.

Now the balance automatically performs its own calibration.

After calibration is completed, the display indicates:

100.00 g in the PE360,
1000.0 g in the PE3600.

Note: if no weight or a weight different from 100 g / 1000 g is placed on the pan while the balance is displaying "- CAL -", first "no CAL" appears, then zero.


In this case, no calibration has taken place.

Note: Certified balances can be checked, but not calibrated.

OPERATION: Selecting the weight unit


The PE360/3600 balances can display weights in grams and even in several nonmetric units. Thus, they can work for instance as ounce balances. The selected weight unit remains stored even in the case of a power interruption but it can be changed at any time.

Note: In certified balances, the weight unit must be selected before certification is carried out.

Selecting the weight unit

- Disconnect power-line cable (14).
- Press the control bar (1) and hold it down. With the bar pressed, connect power-line cable (14).
- The display indicates now "Unit". On the right of the display, the weight units light up in succession.
- When the desired weight unit lights up, release the control bar (1).
- The display indicates all the numerals and characters (in the same way as when being switched on) and then follows the zero display.


The number of decimal places depends on the balance model and the weight unit.

- If during the changing of the weight units a given unit does not appear on the display, that particular unit can not be selected, e.g. the kilogram on the PE360.

Note: Calibration is always operated in the same way, regardless of the selected weight unit. Calibration weight must be 100 g (PE360), or 1000 g (PE3600).

The displayed value can be calculated from the table on the right. Example of the PE360 as an ounce balance: 100 g \triangleq 3.5274 oz.


Conversion factors

Ounce	1 oz	\triangleq	28.349523125 g
	1 g	\triangleq	0.035273962 oz
Pound	1 lb	\triangleq	453.59237 g
	1 g	\triangleq	0.002204623 lb
Pennyweight	1 dwt	\triangleq	1.55517384 g
	1 g	\triangleq	0.643014931 dwt
Troy ounce	1 ozt	\triangleq	31.1034768 g
	1 g	\triangleq	0.032150747 ozt
Grain	1 GN	\triangleq	0.06479891 g
	1 g	\triangleq	15.43235835 GN
Carat	1 ct	\triangleq	0.2 g
	1 g	\triangleq	5 ct
Tael	1 tl	\triangleq	37.4375 g
	1 g	\triangleq	0.026711185 tl

OPERATION: How weight is determined or tared out

The PE360 DeltaRange and the PE3600 DeltaRange allow you to weigh very easily in two different ranges: either the coarse range, or the fine range whose readability is 10 times greater and which can be shifted to any place within the entire balance weighing range. This fine range extends from the zero display in both plus and minus direction. Besides, the PE balance measures one decimal place more than it displays, thus increasing the display accuracy of PE balances even more. Taring means setting the display to zero again while a given weight (e.g. a container) is on the weighing pan. In this way, it is possible to always weigh in starting with the zero display. Only the weight of the weighed-in object appears in the display. Of course, combined container weight and object to be weighed may not exceed the balance weighing range.

Weighing


- Briefly press the control bar (1); the rest the display to zero.
- Place object on pan.
- Read weight display when stability detector lamp (15) no longer lights.

The fine range is available starting from the zero position, i.e. the display indicates an additional decimal place.

If the fine range is exceeded, this decimal place is blanked out: the balance now displays only the coarse range.

Switching back to the fine range:

- Press control bar (1); the balance is also set to zero.

Display accuracy (result rounding)

Your balance always measures one place more than it displays. This last measured digit is rounded off according to the 4/5 principle.

Example PE3600:


The balance measures 94.314 g;
it displays 94.31 g.

or

The balance measures 94.315 g
it displays 94.32 g

Overload indication

If the balance weighing range is exceeded, the display blanks out except for the upper horizontal segments of the digits. This is how the balance indicates overload.

Taring

- Place a container on the pan: its weight is displayed.
- Briefly press the control bar (1), i.e. tare: The balance displays zero.

The balance weighing range – minus the container weight – is now available for weighing.


By taring, the fine range has been shifted to the current zero position

If bar (1) is pressed while there is instability – i.e., while the pilot of the stability detector is on –, the entire display goes out until stability is achieved; then the balance displays zero

OPERATION: How to weigh in

Weighing-in means filling loose or liquid objects into a container on the balance until a desired weight (the target weight) has been reached.


PE balances are equipped with DeltaDisplay which makes it easier to weigh in because, when fast weight changes occur, the last number of the display is temporarily blanked out.

Weighing-in

- Place the container on the weighing pan.
- Tare by pressing control bar (1): the balance displays zero.
- Fill in objects to be weighed up to the desired target weight.

Note: While weight is changing in the display, the stability detector lamp (15) is on. It only goes off when the result is stable.


Reading results

While filling in steadily, pay attention to the first two digits (16) only until in the range of the target weight.

While filling in very slowly towards the target weight, pay carefully attention to the digits to the right (17).

When weight is increasing quickly, the last digit (18) blanks out temporarily; when weighing in slowly, this digit reappears, so far as the fine range has not been exceeded.

Weighing in several components

If several components must be weighed in, one after the other, each weighing can be tared out and the next weighing can proceed again from zero and with the fine range readability, until the combined weight of container and objects reaches the end of the weighing range:


- in the PE360: 360 g
- in the PE3600: 3600 g

OPERATION: How to make subtractive weighings; how to read deviations from a reference weight

A container filled with a substance can be tared out. This makes it possible to weigh out from the container without having to calculate how much substance has been removed.


Deviations from a reference weight can be read with the correct sign in front, provided the reference weight has been tared out beforehand. In any case, the balance display has to be switched on before a load (e.g., the reference weight) is placed on the pan.

The deviation is displayed with the fine range readability so far as the reference weight does not exceed the fine range.

Weighing out


- Place container with the substance on the weighing pan.
- Tare; the balance display zero.
- Remove some of the substance; the weight of the removed substance appears in the display with a negative sign in front.

Setting a reference weight

- Place reference weight on pan (weight or master).
- Tare: the balance displays zero.
- Remove reference weight; the balance now displays the reference weight with a negative sign in front, e.g. - 300.00 g

Now, deviations from this reference weight can be read directly.

Determining a deviation from a reference weight

- Place object to be compared against reference on pan.

If the object is heavier than the reference, the deviation appears in the display with a positive sign in front because:

$$- 300.00 \text{ g} + 300.62 \text{ g} = + 0.62 \text{ g}$$

If the object is lighter than the reference, the deviation appears with a negative sign in front, because:

$$- 300.00 \text{ g} + 299.42 \text{ g} = - 0.58 \text{ g}$$


FURTHER CAPABILITIES OF PE BALANCES: Level indicator Field Installation Kit, glass draft shield, GD hanger

There is a Level indicator Field Installation Kit available for PE balances that are frequently moved from location to location.

This allows the balance to be leveled at every new location. Thus, the accuracy of the balance remains constant.

A glass draft shield is available as an option for the PE360.


A GD hanger is built in as standard equipment in the PE balances; this is used for weighing below the balance.

Level indicator Field Installation Kit


The Field Installation Kit contains two foot screws (28) and one level indicator (27). Installation can be carried out either by the user himself, or by METTLER TOLEDO Service.

Please note: The balance should be leveled before being calibrated. This way, the leveling kit can perform its function.

Glass draft shield (for PE360)

With the PE360, it is possible to replace the plastic draft shield that comes with the balance with the METTLER TOLEDO glass draft shield. Use of this shield prevents the weighing result from being affected when using tall weighing containers.

GD hanger

For weighings below the balance, there is a GD hanger built into the balance. Here, the object to be weighed is hooked onto the weighing cell; this requires an opening in the weighing table.


- Remove weighing pan and pan support; tip balance to its side (do not place balance on it's upper part).
- Remove cover (23).
- Hook hanger for object to be weighed onto hook (24).
- Place balance on its feet; install pan support and weighing pan.
- Press tare with hanger assembled; proceed with weighing.

FURTHER CAPABILITIES OF PE BALANCES: Application Packages, Data Output

A Data Output is available for the PE balances. This can be installed by the user very easily. For special applications, there are accessories available for PE balances; a selection of these is presented on the following two pages.

Your routine weighings can be greatly simplified by using METTLER TOLEDO Application Packages which are programmed for a range of particular weighing operations.

They contain a GE305 Application Input Device or a GE310 Terminal and a Data Output board.

Data Output


The following peripheral instruments can be connected to the data output socket (2):

- GA24 Calculating printer
- GA44 Thermoprinter
- METTLER TOLEDO CL instruments
- Non-METTLER TOLEDO instruments with CL or RS232C data interfacing.

To the socket (3) can be connected:

- GE305 Application Input Device
- GE310 Terminal
- GE48 Auxiliary Display.

The Data Output alone (without Application Package) can be delivered with the balance or ordered at a later date, see "ACCESSORIES".


Application Packages

a) with the GE305

Application Input Device:

- LabPac which offers interesting application functions for laboratories.
- ProPac which offers important application functions for production.
- GemPac I/II for weighing precious metals or stones.


b) with a GE310 Terminal:

- StatPac for statistical quality control (SQC11).
- CountPac for piece counting with the possibility of entering the reference weight, a pretare value etc.
- DataPac for designing any desired application for use with a computer. The GE54 Output Module is available as an option; this instrument makes control of a variety of additional functions possible.

PERIPHERAL INSTRUMENTS: Thermoprinter, Infrared Dryer, Auxiliary Display

The GA44 Thermoprinter, presented below, can be connected to any PE balance if it is equipped with a data output (also see the previous page).

The Thermoprinter, and also the GE48 Auxiliary Display, can be connected to the balance at the same time as the GE305 Application Input Device or the GE310 Terminal.


GA44 Thermoprinter

For hard-copy printout of weighing results.

Maintenance-free thermoprinter (dot printer without ink ribbon); compact, fast and quiet.

Non-fading paper, not sensitive to print or pressure; rolls 30 meters long for approx. 5000 lines of print; tape width: 57 mm (2 1/4").

Separate power supply: 110/130/220/240 V
Dimensions (W x D x H): 185x200x60 mm


LP15B Infrared Dryer

For determination of moisture content or of dry quantity in samples.

If the GE305 is used along with this unit, the balance will display results in either % or g, as desired.

Separate power supply:
110/220 V, 140 VA.
Dimensions (W x D x H): 190x200x230 mm
12 temperature steps
(approx. 50...160 °C).








GE48 Auxiliary Display

To make the work station more efficient. The Auxiliary Display can be attached to the balance, or it can stand alone. Power supply from the balance. Balance Data Output is required.





ACCESSORIES

Optional equipment:

	Order No.
 Glass draft shield (for PE360 only)	43970
 Calibration weight 100 g for PE360 1000 g for PE3600	47900 47905
 Data Output, 016 Option with RS232C cable set Data Interface, 017 Option with CL cable set	48330 59817
 CL cable (5 poles) 1.5 m GA44 cable (7 poles) 1.5 m RS232C cable 1.5 m Adapter for foot pedal / hand key	47936 47926 47927 47473
 Level indicator kit	47659
 Microfuses (set of 3) 125 mA slow blowing (95...120 V) 63 mA slow blowing (190...240 V)	26172 46328
 Animal weighing pan (for PE3600 only)	47644
 Application Packages: for labs for production for fillings operations (SQC11) for counting for gem weighing (PE360) for gem weighing (PE3600) for applications with computer	LabPac ProPac StatPac CountPac GemPac I GemPac II DataPac

 Auxiliary display Overhead display	GE48 GE53
 Thermoprinter Calculating printer	GA44 GA24
 Infrared dryer 110/220 V, 140 VA	LP15B


Standard equipment:

	Order No.
 Screwdriver	50279
 In-use cover	47646
 Retainer ring (with PE3600)	41156
 Weighing pan	PE360 47085 PE3600 47087
 Pan support	PE360 47081 PE3600 47083
 Power cable	Neutral 87576 Switzerland 87920 Germany 87925 USA 88668
 Plastic draft shield with 2 sections (for PE360 only)	47910

MAINTENANCE: Installing the in-use cover and changing the fuse

To protect the balance against dirt, damage and caustic chemical substances, the balance is shipped with an installed in-use cover. The in-use cover is made of hard plastic and cannot be folded. Operation of the balance is in no way hindered by the in-use cover.


In case the microfuse needs to be replaced, there is a spare fuse in the fuse holder. If the microfuse needs to be replaced repeatedly for no obvious reason, contact METTLER TOLEDO Service.

Removing the in-use cover


- Remove weighing pan and pan support.
- Turn retaining ring (25) (in the PE360, the draft shield base (19)) counterclockwise until it snaps free and can be lifted out.
- Remove in-use cover (26).
- Place pan support and weighing pan back on balance.

(The retaining ring is installed in the PE3600 only if an in-use cover is used).

Installing the in-use cover

- Remove weighing pan and pan support.
- First hook up in-use cover on the rear of balance housing and then pull it tight at the balance front.
- Place retaining ring (25) (in the PE360, the draft shield base) then turn clockwise to stop position.
- Place pan support and weighing pan back on balance.


Replacing the microfuse

- Disconnect the power cable.
- Remove the fuse holder (4) with a screwdriver.
- Remove defective fuse and install a new fuse (Spare fuse is in the fuse holder).
Fuse values:
125 mA slow-blowing with 95...120 V
63 mA slow-blowing with 190...240 V
- Reinsert fuse holder.
- Plug power cable back in.

MAINTENANCE: Cleaning and calibration

Your balance requires only very little maintenance: it should be cleaned regularly and calibrated periodically. Before checking calibration, the balance should be left connected to the power source for 30 minutes (warm-up time). If your balance is equipped with a leveling kit, it should be leveled before it is calibrated.


Cleaning

Please note: The balance must not be placed on its top, because that would damage the measuring cell.


As needed, the weighing pan, in-use cover, the draft shield (PE360), and if necessary, the housing should be cleaned.

If possible, do not use any strong solvents, so as not to damage the paint of the balance housing.

If you own both PE360 and PE3600 balances, please make sure not to interchange the pan supports and weighing pans of the various balances after cleaning.

The pan support of the PE360 (soft springs) must not be used on a PE3600 (hard springs) or the other way around.

Calibration

To have your balance always provide accurate weighing results, it must be correctly calibrated.

Check calibration regularly, e.g. with continuous operation, once a month.

The test weight used with your balance model is listed under "ACCESSORIES, optional".

WHAT'S WRONG IF ...

- ... the entire display does not light up?
- ... only the upper display segments light up?
- ... only the lower display segments light up?
- ... the display indicates "OFF"?
- ... the weighing result is instable?
- ... the weighing result is obviously incorrect?
- ... "no CAL" lights up when calibrating?
- ... the balance displays unclear characters or the balance is blocked?
- ... the balance displays "ERROR"?

Then ...

- The display has not been switched on. Press control bar.
- The power cable has not been plugged in. No power reaching balance.
- The microfuse is defective (see "MAINTENANCE: Changing the microfuse").
If this happens repeatedly: check the voltage setting and the value of the microfuse.
If both are correct, contact METTLER TOLEDO Service.
- Too heavy weight has been placed on the pan (weighing range has been exceeded).
- Balance is defective (contact METTLER TOLEDO Service).
- Pan support or weighing pan not installed.
- In-use cover touching the pan support (retaining ring for in-use cover not installed).
- Balance switched on before unloading the pan. Tip: tare with no load on pan.
- Power supply was temporarily interrupted: check both power cable connectors for tight fit.
Press control bar to switch on balance again.
- Slingle control bar was not pressed down completely when balance was switched on.
- Drafts (PE360: use draft shield).
- The weighing table is instable.
- The object on pan is instable, e.g., live animal.
- Jumper plug on wrong voltage selection position (see "PRELIMINARY STEPS: checking the voltage setting").
- The balance has been placed on an uneven surface.
- Balance is out of calibration (see "OPERATION: How to calibrate the balance").
- The in-use cover is touching the pan support (install the retaining ring or draft shield base).
- Tare was not pressed before weighing was made (control bar not pressed).
- Object on pan touching draft shield (if installed).
- An incorrect pan support has been installed on the PE3600.
- An improper calibration weight has been used (for weights, see page 7).
- The electronic system is not working properly. Unplug power cable, then plug it in again.
- There is a defect in the balance electronics; contact METTLER TOLEDO Service.

TECHNICAL DATA

	PE360	PE3600
Weighing range	360 g	3600 g
- in METTLER TOLEDO DeltaRange	60 g	600 g
Readability	0.01 g	0.1 g
- in METTLER TOLEDO DeltaRange	0.001 g	0.01 g
Tare range (subtractive)	360 g	3600 g
Admissible ambient temperature (during operation)		
- Temperature	0...+ 40 °C	
- Altitude	-500...+6000 m	
- Humidity (non-condensing)	15...85 %	
- Vibrations	0.3 m/s ²	
Reproducibility (s)	0.003 g	0.03 g
- in METTLER TOLEDO DeltaRange	0.001 g	0.01 g
Linearity	± 0.01 g	± 0.1 g
- in METTLER TOLEDO DeltaRange	± 0.002 g	± 0.02 g
Stabilization time (typical)	2.5 s	2.5 s
Display change	0.2 s	0.2 s
Sensitivity drift (10...30 °C)	± 4•10 ⁻⁶ /°C	± 4•10 ⁻⁶ /°C
Result deviation (with balance inclined by 1:1000)	± 0.003 g	± 0.03 g
Power supply:		
- Voltage adjustable	100 V / 115 V / 200 V / 230 V	
- Tolerance	+ 10 % / - 15 %	
- Frequency	50...60 Hz	
- Power consumption	approx. 6 VA	
Weighing pan (chrome-nickel steel)	Ø 130 mm	Ø 170 mm
Housing dimensions (W x D x H)	195 x 310 x 67 mm	195 x 310 x 67 mm
Weight	3.7 kg	3.7 kg

Data Output (016 Option)

Combined current loop and RS232C data interface. 2-wire data interface: unidirectional, galvanically separated. Characters: alphanumeric; character set: ASCII/ISO, 7-bit even parity. Transfer speed: 300/2400 baud. Plug connector: 15-pole, type MiniMettler. CL mode: (for METTLER TOLEDO CL instruments and non METTLER TOLEDO instruments with active current loop input). Current loop: 20 mA (passive). RS232 mode: (for non-METTLER TOLEDO instrument) voltage-controlled data interface conforming to the RS232C/V24 standard.

Data Interface (017 option)

Current loop data interface for METTLER TOLEDO CL instruments and non-METTLER TOLEDO instrument with active current loop input. Loop current 20 mA (passive). Four-wire data interface: bidirectional, galvanically separated. Character alphanumeric. Character code: ASCII/ISO, 7-bit, selectable parity. Transfer: asynchronous with and without handshake. Transfer speed: 300/600/1200/2400/4800/9600 baud. Connector socket: 15-pole.

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© Mettler-Toledo GmbH 2001 701344 Printed in Switzerland 0112/2.33
Mettler-Toledo GmbH, Laboratory & Weighing Technologies, CH-8606 Greifensee, Switzerland